

Study Guide for Lab Exam 2017

Basic Lab Report Requirements:

NO PERSONAL PRONOUNS

-Materials should be organized (list or columns)

-Number the steps of the procedure

-Be specific in the procedure

-Procedure is written in present tense and starts with a verb

-Record results and repeat steps (where needed)

NO PERSONAL PRONOUNS

→ Data needs to be presented in an organized, neat manner (use a table with a title)

-Restate your aim and summarize test results in the conclusion

NO PERSONAL PRONOUNS

*2 sources of error (that are not silly mistakes: mislabelling or miscalculating)

Lab Tests:

You should **KNOW** the information in the table as it will not be provided to you in the final (we have done these tests several times)

| Test | Method | Used For | Possible Results |
|-----------------------|---|--|--|
| Conductivity | -The electrodes of a conductivity meter are placed into 2 drops of the unknown | -Detecting electrolytes, <u>minerals</u> and/or <u>salts</u> in a liquids | -Light flashes quickly -Light flashes slowly -No light |
| Cobalt Chloride (CCP) | -A piece of cobalt chloride paper is placed in one drop of the unknown -Record any change in colour the is <u>immediate</u> | -Detecting the presence of <u>water</u> -Change in colour means water is present | -Paper <u>turns pink</u> -Paper stays dark blue |
| pH (RLP and BLP) | - Blue litmus paper is placed in the well of plate and one drop of the unknown is dropped on top: acids are detected - Red litmus paper is placed in the well of plate and one drop of the unknown is dropped on top: bases are detected | -If the blue paper turns red it means it is acidic -If the red paper turns blue it means it is basic -If neither changes, the pH is close to 6-7, which is neutral -Both paper tests need to be done in order to determine the pH | -BLP - <u>Turns red</u> - Stays blue -RLP - <u>Turns blue</u> - Stays red |
| Protein Presence | -Make a control by placing 3 drops of each NaCO_3 and CuSO_4 -Then repeat this and add 6 drops of the unknown -Wait at least <u>2 min</u> to observe results | -Detecting the presence of protein -Proteinuria -Protein in blood/urine | -Stays blue - <u>Turns violet/dark blue/green/ around the edges</u> (noticeable color change must occur) |

Procedure

Results

Place
Put
Record } action word
headings use tests in title

-ruler

* Make a control *

* if prote

| | | | |
|-----------------------|--|--|--|
| Glucose Presence * | -A glucose strip is placed in the in the well with tweezers and one drop of the unknown is added - <u>wait 30 seconds</u> for observations | -Detecting the presence of carbohydrates (colour change) <u>Glucose/Sugar</u> | - <u>Turns green</u> (glucose) -Stays yellow (o glucose) |
| Salt presence | - <u>Conductivity</u> meter electrodes are placed in the unknown | -Detecting presences of electrolytes | - <u>More than 2 lights</u> will turn on from the conductivity meter (salt) -0 or 1 light goes on |
| Starch presence | -Lugol's Solution -place 5 drops of Unknown A to a clean, dry well of the spot plate -Add 2 drops of Lugol's solution to the well containing the unknown using the pipette -Stir the solution with a clean stirring rod -Observe if unknown turns <u>green/black</u> and <u>record</u> observations. | -Detecting the presence of starch and determining if amylase (saliva) will have an effect on the substance | - <u>Unknown turns</u> green/black/brown (starch) -Turned <u>yellow/orange</u> (no starch) |
| Fat Presence | -Add 20 drops of Unknown A to the "A" test tube -Using dropper bottle, add 6 drops of ethanol to the "A" test tube -Place thumb over mouth of test tube and shake vigorously side-to-side -Observe for an emulsion (<u>white cloudy</u> appearance) | - Emulsion with Ethanol | - <u>Unknown turns</u> clumps, is cloudy, emulsifies (fat) - <u>Stays clear</u> (No fat) - <u>Stays cloudy</u> |